
Product Information

Product SKU:	CAB19403	Gene ID:	9804	Size:	20uL, 100uL
Clone No:	ARC5002-01	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	16kDa/	Conjugate:	Unconjugated
Calculated MW:	16kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein-transporting ATPase activity and unfolded protein binding activity. Involved in protein targeting to mitochondrion. Located in mitochondria-associated endoplasmic reticulum membrane and mitochondrial outer membrane.

Recommended Dilution: WB,1:2000 - 1:6000 IHC-P,1:50 - 1:200 IF/ICC,1:100 - 1:500 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells

Synonyms: MAS20; MOM19; TOM20

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOM20 (NP_055580.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.