

CAB19421

Product Information

Product SKU:	CAB19421	Gene ID:	4839	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	89kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables RNA binding activity. Involved in positive regulation of cell population proliferation; regulation of signal transduction by p53 class mediator; and ribosomal large subunit assembly. Located in nucleolus.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	NOL1; p120; NSUN1; NOP120; NOP2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 355-590 of human NOP2 (NP_001245237.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.