

PTMA Rabbit Polyclonal Antibody



CAB1956

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

12kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5757

Uniprot

P06454

Synonyms:

PTMA; TMSA

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human PTMA (NP_002814.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
