

CAB19587

Product Information

Product SKU:	CAB19587	Gene ID:	8887	Size:	20uL, 100uL
Clone No:	ARC2184	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	91kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: T6BP; TXBP151; CALCOCO3; TRAF6BP/TAX1BP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 690-789 of human TRAF6BP/TAX1BP1 (Q86VP1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.