

PAK1 Rabbit Monoclonal Antibody



CAB19608

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

66KDa

Calculated MW:

66kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a family member of serine/threonine p21-activating kinases, known as PAK proteins. These proteins are critical effectors that link RhoGTPases to cytoskeleton reorganization and nuclear signaling, and they serve as targets for the small GTP binding proteins Cdc42 and Rac. This specific family member regulates cell motility and morphology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2010]

Immunogen information

Gene ID:

5058

Uniprot

Q13153

Synonyms:

PAK1; PAKalpha

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

A synthesized peptide derived from human PAK1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images

Western blot - PAK1 Rabbit mAb (CAB19608)

