

Secretagogin (SCGN) Rabbit Monoclonal Antibody

CAB19615



Product Information

Product SKU:	CAB19615	Gene ID:	10590	Size:	20uL, 100uL
Clone No:	ARC2196	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	SEGN; CALBL; SECRET; setagin; DJ501N12.8; Secretagogin (SCGN)
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 177-276 of human Secretagogin (Secretagogin (SCGN)) (O76038).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.