

CAB19628

Product Information

Product SKU:	CAB19628	Gene ID:	27	Size:	20uL, 100uL
Clone No:	ARC2200	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	128kDa	Conjugate:	Unconjugated
Calculated MW:	128kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the Abelson family of nonreceptor tyrosine protein kinases. The protein is highly similar to the c-abl oncogene 1 protein, including the tyrosine kinase, SH2 and SH3 domains, and it plays a role in cytoskeletal rearrangements through its C-terminal F-actin- and microtubule-binding sequences. This gene is expressed in both normal and tumor cells, and is involved in translocation with the ets variant 6 gene in leukemia. Multiple alternatively spliced transcript variants encoding different protein isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ARG; ABLL; ABL2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human ABL2 (P42684).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.