

Product Information

Product SKU:	CAB19642	Gene ID:	957	Size:	20uL, 100uL
Clone No:	ARC2210	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	48kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is similar to E-type nucleotidases (NTPases)/ecto-ATPase/apyrases. NTPases, such as CD39, mediate catabolism of extracellular nucleotides. ENTPD5 contains 4 apyrase-conserved regions which is characteristic of NTPases.

Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	PCPH; CD39L4; NTPDase-5; ENTPD5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human ENTPD5 (O75356).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.