

## CAB19693

---

### Product Information

<b>Product SKU:</b>	CAB19693	<b>Gene ID:</b>	596	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC0173	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	26kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	26kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IP,0.5µg-4µg antibody for 400µg-600µg extracts of whole cells

**Synonyms:** Bcl-2; PPP1R50

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human Bcl-2 (P10415).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.