

CAB19727

Product Information

Product SKU:	CAB19727	Gene ID:	4250	Size:	20uL, 100uL
Clone No:	ARC2255	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	10kDa	Conjugate:	Unconjugated
Calculated MW:	10kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in androgen receptor signaling pathway. Predicted to be located in extracellular region. Predicted to be active in extracellular space.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	MGB1; UGB2; PSBP1; Mammaglobin A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human Mammaglobin A (Q13296).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.