

---

## Product Information

<b>Product SKU:</b>	CAB19751	<b>Gene ID:</b>	5962	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2277	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	80kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	69kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** DFNB24; Radixin

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Radixin (P35241).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.