
Product Information

Product SKU:	CAB19809	Gene ID:	84864	Size:	20uL, 100uL
Clone No:	ARC2336	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	53kDa	Conjugate:	Unconjugated
Calculated MW:	53kDa	Isotype:	IgG

Immunogen Information

Background: MINA is a c-Myc (MYC; MIM 190080) target gene that may play a role in cell proliferation or regulation of cell growth. (Tsuneoka et al., 2002 [PubMed 12091391]; Zhang et al., 2005 [PubMed 15897898]).

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: ROX; MDIG; MINA; NO52; JMJD10; MINA53; RIOX2/MINA53

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human RIOX2/MINA53 (Q8IUJ8).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.