

CAB19817

Product Information

Product SKU:	CAB19817	Gene ID:	145258	Size:	20uL, 100uL
Clone No:	ARC2343	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SAMS; GSC
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GSC (P56915).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.