

CAB19856

Product Information

Product SKU:	CAB19856	Gene ID:	29982	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Involved in autophagy. Located in cytoplasm. Colocalizes with phosphatidylinositol 3-kinase complex, class III.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: COPR; COPR1; COPR2; NRBF-2; F2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 10-220 of human NRBF2 (NP_110386.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.