

[KO Validated] IPO9 Rabbit Polyclonal Antibody



CAB19860

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

130KDa

Calculated MW:

115kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

55705

Uniprot

Q96P70

Synonyms:

IPO9; Imp9

Immunogen:

Recombinant protein of human IPO9.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and IPO9 knockout (KO) HeLa cells, using IPO9 antibody (CAB19860) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.