

**[KO Validated] PSME4 Rabbit Polyclonal
Antibody**
CAB19873



Product Information
Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

211KDa

Calculated MW:

24kDa/113kDa/139kDa/211kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

23198

Uniprot

Q14997

Synonyms:

PSME4; PA200

Immunogen:

Recombinant protein of human PSME4.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and PSME4 knockout (KO) HeLa cells, using PSME4 antibody (CAB19873) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.