

## CAB19883

---

### Product Information

<b>Product SKU:</b>	CAB19883	<b>Gene ID:</b>	7410	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	95kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	101kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** VAV-2; V2

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 570-700 of human VAV2 (NP\_001127870.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.