

CAB19888

Product Information

Product SKU:	CAB19888	Gene ID:	10147	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	145kDa	Conjugate:	Unconjugated
Calculated MW:	120kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the arginine/serine-rich family of splicing factors. The encoded protein functions in mRNA processing. Alternatively spliced transcript variants have been described.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	SFRS14; SRFS14; P2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human SUGP2 (NP_055699.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.