

[KO Validated] SUGP2 Rabbit Polyclonal Antibody

CAB19888



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

145KDa

Calculated MW:

115kDa/116kDa/120kDa

Applications:

WB

Reactivity:

Human, Mouse

Protein Background

This gene encodes a member of the arginine/serine-rich family of splicing factors. The encoded protein functions in mRNA processing. Alternatively spliced transcript variants have been described.

Immunogen information

Gene ID:

10147

Uniprot

Q8IX01

Synonyms:

SUGP2; SFRS14

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant protein of human SUGP2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and SUGP2 knockout (KO) HeLa cells, using SUGP2 antibody (CAB19888) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.