

**[KO Validated] TXLNA Rabbit Polyclonal
Antibody
CAB19889**



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

80KDa

Calculated MW:

61kDa

Applications:

WB

Reactivity:

Human

Protein Background

Immunogen information

Gene ID:

200081

Uniprot

P40222

Synonyms:

IL14; TXLN; TXLNA

Immunogen:

Recombinant protein of human TXLNA.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and TXLNA knockout (KO) HeLa cells, using TXLNA antibody (CAB19889) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.