

**[KO Validated] NSDHL Rabbit Polyclonal
Antibody
CAB19896**



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

38KDa

Calculated MW:

41kDa

Applications:

WB IF

Reactivity:

Human, Rat

Protein Background

The protein encoded by this gene is localized in the endoplasmic reticulum and is involved in cholesterol biosynthesis. Mutations in this gene are associated with CHILD syndrome, which is a X-linked dominant disorder of lipid metabolism with disturbed cholesterol biosynthesis, and typically lethal in males. Alternatively spliced transcript variants with differing 5' UTR have been found for this gene.

Immunogen information

Gene ID:

50814

Uniprot

Q15738

Synonyms:

NSDHL; H105E3; SDR31E1; XAP104

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human NSDHL.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
