

CAB19901

Product Information

Product SKU: CAB19901	Gene ID: 8030	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 68kDa	Conjugate: Unconjugated
Calculated MW: 53kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a coiled-coil domain-containing protein. The encoded protein is ubiquitously expressed and may function as a tumor suppressor. A chromosomal rearrangement resulting in the expression of a fusion gene containing a portion of this gene and the intracellular kinase-encoding domain of the ret proto-oncogene is the cause of thyroid papillary carcinoma.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: H4; PTC; TPC; PTC1; TST1; D10S170; C6

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 55-222 of human CCDC6 (NP_005427.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.