

# [KO Validated] HAGH Rabbit Polyclonal Antibody



CAB19915

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

30KDa

### Calculated MW:

25kDa/28kDa/33kDa

### Applications:

WB IHC IF

### Reactivity:

Human, Mouse, Rat

## Protein Background

The enzyme encoded by this gene is classified as a thiolesterase and is responsible for the hydrolysis of S-lactoyl-glutathione to reduced glutathione and D-lactate. Three transcript variants encoding different isoforms have been found for this gene.

## Immunogen information

### Gene ID:

3029

### Uniprot

Q16775

### Synonyms:

HAGH; GLO2; GLX2; GLXII; HAGH1

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50  
- 1:200 IF 1:10 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:


Recombinant protein of human HAGH.

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of extracts from normal (control) and HAGH knockout (KO) HeLa cells, using HAGH antibody (CAB19915) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.