

CAB19920

Product Information

Product SKU:	CAB19920	Gene ID:	22929	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	41kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes an enzyme that synthesizes selenophosphate from selenide and ATP. Selenophosphate is the selenium donor used to synthesize selenocysteine, which is co-translationally incorporated into selenoproteins at in-frame UGA codons.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SPS; SELD; SPS1; S1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human SEPHS1 (NP_036379.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.