

**[KO Validated] RAB1B Rabbit Polyclonal
Antibody**
CAB19935



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

22KDa

Calculated MW:

22kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

81876

Uniprot

Q9H0U4

Synonyms:

RAB1B

Immunogen:

Recombinant protein of human RAB1B.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02%
sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and RAB1B knockout (KO) HeLa cells, using RAB1B antibody (CAB19935) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.