

CAB19944

Product Information

Product SKU:	CAB19944	Gene ID:	166	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is similar in sequence to the amino terminus of Drosophila enhancer of split groucho, a protein involved in neurogenesis during embryonic development. The encoded protein, which belongs to the groucho/TLE family of proteins, can function as a homooligomer or as a heterooligomer with other family members to dominantly repress the expression of other family member genes. Three transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: AES; GRG; ESP1; GRG5; AES-1; AES-2; Grg-5; ES

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human AES (NP_945320.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.