

[KO Validated] STRAP Rabbit Polyclonal Antibody



CAB19989

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

39KDa

Calculated MW:

38kDa/39kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

11171

Uniprot

Q9Y3F4

Synonyms:

STRAP; MAWD; PT-WD; UNRIP


Immunogen:

Recombinant protein of human STRAP.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and STRAP knockout (KO) 293T cells, using STRAP antibody (CAB19989) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1s.