

**CAB19989**

---

**Product Information**

<b>Product SKU:</b>	CAB19989	<b>Gene ID:</b>	11171	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	39kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	38kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** Enables RNA binding activity. Involved in maintenance of gastrointestinal epithelium; negative regulation of transforming growth factor beta receptor signaling pathway; and spliceosomal snRNP assembly. Located in cytosol. Part of SMN complex. Implicated in adenocarcinoma; colorectal carcinoma; large cell carcinoma; lung carcinoma; and squamous cell neoplasm.

**Recommended Dilution:** WB,1:1000 - 1:2000

**Synonyms:** MAWD; PT-WD; UNRIP; AP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human STRAP (NP\_009109.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.