

**Os03g0203200/IOC\_Os03g10620 Polyclonal  
Antibody  
CAB20012**



---

**Product Information**

**Size:**

50uL, 100uL, 200uL

**Observed MW:**

**Calculated MW:**

**Applications:**

WB

**Reactivity:**

Oryza sativa

**Antibody Information**

**Recommended dilutions:**

WB 1:500 - 1:2000

**Source:**

Rabbit

**Isotype:**

IgG

**Purification:**

Affinity purification

**Protein Background**

**Immunogen information**

**Gene ID:**

4331983

**Uniprot**

P42224

**Synonyms:**

D14; D88; HTD2

**Immunogen:**

Recombinant protein of Oryza sativa  
Os03g0203200/IOC\_Os03g10620

**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---