

POSTN Rabbit Polyclonal Antibody



CAB2003

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

80kDa/83kDa/87kDa/90kDa/
93kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms.

Immunogen information

Gene ID:

10631

Uniprot

Q15063

Synonyms:

POSTN; OSF-2; OSF2; PDLPOSTN; PN; periostin

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 472-751 of human POSTN (NP_001129408.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
