

## CAB20139

---

### Product Information

<b>Product SKU:</b>	CAB20139	<b>Gene ID:</b>	81558	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	48kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	48kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** FAM117A (Family With Sequence Similarity 117 Member A) is a Protein Coding gene. An important paralog of this gene is FAM117B.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 110-240 of human FAM117A (NP\_110429.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.