

TMEM106B Rabbit Polyclonal Antibody

CAB20165



Product Information

Product SKU:	CAB20165	Gene ID:	54664	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	31kDa	Isotype:	IgG

Immunogen Information

Background: Enables ATPase binding activity. Involved in dendrite morphogenesis and lysosome localization. Located in endosome and lysosomal membrane. Implicated in hypomyelinating leukodystrophy.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	HLD16; TMEM106B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 148-274 of human TMEM106B (NP_060844.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.