

CAB20169

Product Information

Product SKU:	CAB20169	Gene ID:	55355	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	80kDa	Conjugate:	Unconjugated
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background:	Enables histone binding activity and identical protein binding activity. Involved in several processes, including CENP-A containing chromatin assembly; regulation of DNA binding activity; and regulation of protein-containing complex assembly. Located in kinetochore; mitochondrion; and nuclear lumen.
Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	FAKTS; URLC9; hFLEG1; HJURP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human HJURP (NP_060880.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.