

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB20255 | Gene ID: | 23246 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|--------|------------|--------------|
| Observed MW: | 110kDa | Conjugate: | Unconjugated |
| Calculated MW: | 84kDa | Isotype: | IgG |

Immunogen Information

Background: Enables RNA binding activity. Involved in regulation of cell cycle; regulation of signal transduction by p53 class mediator; and ribosomal large subunit biogenesis. Located in chromosome; nucleolus; and nucleoplasm. Part of PeBoW complex.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 600-680 of human BOP1 (NP_056016.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.