

Product Information

Product SKU:	CAB20310	Gene ID:	80216	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	108kDa/130kDa	Conjugate:	Unconjugated
Calculated MW:	139kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene encodes an alpha kinase. Mice which were homozygous for disrupted copies of this gene exhibited coordination defects (PMID: 21208416). Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	LAK; ROSAH; 8430410J10Rik; ALPK1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 340-500 of human ALPK1 (NP_001095876.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.