

**METTL9 Polyclonal Antibody**

**CAB20313**



**Product Information**

<b>Product SKU:</b>	CAB20313	<b>Gene ID:</b>	51108	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

**Additional Information**

<b>Observed MW:</b>	40kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	37kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

**Background:** Enables protein-L-histidine N-pros-methyltransferase activity. Predicted to be involved in methylation. Is active in endoplasmic reticulum.

<b>Recommended Dilution:</b>	WB,1:500 - 1:1000
<b>Synonyms:</b>	DREV; PAP1; DREV1; CGI-81; hMETTL9; METTL9
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant protein of human METTL9.
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.