

Product Information

Product SKU:	CAB20390	Gene ID:	493869	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	24kDa	Conjugate:	Unconjugated
Calculated MW:	24kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables peroxidase activity. Predicted to be involved in cellular response to oxidative stress. Predicted to be located in endoplasmic reticulum lumen.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	GPx-8; UNQ847; EPLA847; GSHPx-8; GPX8
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 41-209 of human GPX8 (NP_001008398.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.