

CAB20493

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| Product SKU: | CAB20493 | Gene ID: | 153090 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Rat |

Additional Information

| | | | |
|-----------------------|------------------|-------------------|--------------|
| Observed MW: | Refer to Figures | Conjugate: | Unconjugated |
| Calculated MW: | 132kDa | Isotype: | IgG |

Immunogen Information

Background: DAB2IP is a Ras (MIM 190020) GTPase-activating protein (GAP) that acts as a tumor suppressor. The DAB2IP gene is inactivated by methylation in prostate and breast cancers (Yano et al., 2005 [PubMed 15386433]).

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: AIP1; AIP-1; AF9Q34; DIP1/2; DAB2IP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 205-320 of human DAB2IP (NP_619723.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.