

Product Information

Product SKU:	CAB20531	Gene ID:	85463	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	99kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable endoribonuclease activity and mRNA binding activity. Predicted to be involved in RNA phosphodiester bond hydrolysis, endonucleolytic. Predicted to be active in cytoplasmic ribonucleoprotein granule and nucleus.

Recommended Dilution:	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
Synonyms:	MCPIP3; ZC3H12C
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 735-883 of human ZC3H12C (NP_203748.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.