

---

## Product Information

<b>Product SKU:</b>	CAB20557	<b>Gene ID:</b>	23378	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	58kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	51kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables methylated histone binding activity. Involved in several processes, including cellular response to glucose starvation; intrinsic apoptotic signaling pathway by p53 class mediator; and regulation of gene expression. Located in several cellular components, including cytosol; nuclear lumen; and rDNA heterochromatin. Part of chromatin silencing complex.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** NML; KIAA0409; RRP8

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RRP8 (NP\_056139.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.