

---

## Product Information

|                     |          |                      |        |                    |                 |
|---------------------|----------|----------------------|--------|--------------------|-----------------|
| <b>Product SKU:</b> | CAB20567 | <b>Gene ID:</b>      | 56917  | <b>Size:</b>       | 20uL, 100uL     |
| <b>Clone No:</b>    | -        | <b>Host Species:</b> | Rabbit | <b>Reactivity:</b> | Human,Mouse,Rat |

---

## Additional Information

|                       |       |                   |              |
|-----------------------|-------|-------------------|--------------|
| <b>Observed MW:</b>   | 41kDa | <b>Conjugate:</b> | Unconjugated |
| <b>Calculated MW:</b> | 41kDa | <b>Isotype:</b>   | IgG          |

---

## Immunogen Information

**Background:** This gene encodes a homeobox protein and probable transcriptional regulator. The orthologous protein in mouse controls expression of 3-phosphoinositide dependent protein kinase 1, which promotes survival of pancreatic beta-cells.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** MRG2; MEIS3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 325-375 of human MEIS3 (NP\_001287988.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.