

CAB20631

Product Information

Product SKU:	CAB20631	Gene ID:	4340567	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	88kDa	Conjugate:	Unconjugated
Calculated MW:	88kDa	Isotype:	IgG

Immunogen Information

Background: This protein is involved in the pathway starch biosynthesis, which is part of Glycan biosynthesis.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: ALK; SSS2A; SSII-3; OsSSIIa

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsSSIIa.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.