

## CAB20632

---

### Product Information

<b>Product SKU:</b>	CAB20632	<b>Gene ID:</b>	4328075	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	80kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	73kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This protein is involved in the pathway starch biosynthesis, which is part of Glycan biosynthesis.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** BIP; OsBiP1; OsJ\_05115; OJ1442\_E05.10; OsBip1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa OsBip1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.