

Product Information

|              |          |               |         |             |              |
|--------------|----------|---------------|---------|-------------|--------------|
| Product SKU: | CAB20651 | Gene ID:      | 4347354 | Size:       | 20uL, 100uL  |
| Clone No:    | -        | Host Species: | Rabbit  | Reactivity: | Oryza sativa |

Additional Information

|                |                  |            |              |
|----------------|------------------|------------|--------------|
| Observed MW:   | Refer to figures | Conjugate: | Unconjugated |
| Calculated MW: | 20kDa            | Isotype:   | IgG          |

Immunogen Information

**Background:** Plays an important role in cold adaptation at the seedling stage.

|                       |                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------|
| Recommended Dilution: | WB, 1:500 - 1:2000                                                                                  |
| Synonyms:             | OsZIP73; OsTFX1                                                                                     |
| Purification Method:  | Affinity purification                                                                               |
| Immunogen:            | Recombinant protein of Oryza sativa OsZIP73.                                                        |
| Storage:              | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3. |