

Product Information

Product SKU:	CAB20656	Gene ID:	4340408	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	26kDa	Conjugate:	Unconjugated
Calculated MW:	24kDa	Isotype:	IgG

Immunogen Information

**Background:** OsHAL3, a novel light-responsive protein, plays an essential role in photoperiodic control of flowering time in rice, which is probably mediated by forming a complex with Hd1

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	HAL3; PPCDC; OsHAL3; OsJ_20467; OJ1147_D11.5-1
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsHAL3.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.