

Product Information

Product SKU:	CAB20657	Gene ID:	9266154	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	18kDa	Isotype:	IgG

Immunogen Information

Background: Overexpression of OsFD2, a member of the Poaceae FD2 group, produced smaller leaves with shorter plastochrons, suggesting that OsFD2 controls leaf development.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	OsFD2; OsbZIP55
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsFD2.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.