

Product Information

Product SKU:	CAB20660	Gene ID:	4330838	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Transcriptional activator that mediates abscisic acid (ABA) signaling

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	ABF2; AREB1; BZIP23; OsABF2; OsAREB1; OsBZIP23; OJ1004_A11.20; OsbZIP23
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsbZIP23.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.