

Product Information

Product SKU:	CAB20662	Gene ID:	4342214	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	30kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: an ortholog of Arabidopsis AP1 gene

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	RAP1A; MADS15; FDRMADS3; OsMADS15; RMADS215; OJ1567_G09.109-1
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa OsMADS15.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.