

SD1 Polyclonal Antibody

CAB20664



Product Information

Product SKU:	CAB20664	Gene ID:	4325003	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	43kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Key oxidase enzyme in the biosynthesis of gibberellin

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	sd1; GA20; Sd-1; 20ox2; C20ox2; SD1-2E; GA20OX2; Os20ox2; osGA20ox2; SD1
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Oryza sativa SD1.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.