
Product Information

Product SKU:	CAB20670	Gene ID:	4324691	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	120kDa	Conjugate:	Unconjugated
Calculated MW:	120kDa	Isotype:	IgG

Immunogen Information

Background: Receptor kinase involved brassinosteroid (BR) signal transduction. Regulates, in response to BR binding, a signaling cascade involved in plant development, promotion of cell elongation and flowering (Probable). Activates BR signaling by targeting and phosphorylating BSK3, a positive regulator of BR signaling

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: D61; BRI1; OsBRI1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa D61.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.